

Below Knee Shrinker

Application:

1. Using both hands, scrunch down the shrinker to the seam on the bottom and stretch the opening. The seam of the shrinker should be horizontal or from side to side.
2. Pull the shrinker up on the limb and pull up the sides as you travel up the limb.
3. The shrinker should completely cover the knee and go onto the thigh. The shrinker is very snug. If your limb is still healing, it may be difficult to put on initially. It may be easier to put on with assistance from another person (See Figure 1) or with the use of a donning can (4) Make sure the shrinker is pulled taut with no wrinkles. The shrinker must be pulled up all of the way at all times. This may require you to pull up the shrinker several times throughout a day to keep it snug. The shrinker may fall off while you are sleeping. If this happens, just remember to pull it back on as soon as you can.
4. Wear the shrinker at all times, day and night, except when you are bathing yourself or cleaning the shrinker. When you are cleaning the shrinker, a 2nd shrinker or an elastic bandage should be applied to the limb to prevent swelling.
5. Remember to keep your leg straight when sitting or laying down. Do not sit with your knee bent.

Cleaning and Maintenance:

Follow the manufacturer's directions to clean the shrinker. Most recommend hand or machine wash on low settings with mild soap. Machine dry on low or hang to dry.